

Stent Roller

With the Stent Roller, stents in different sizes can be inspected or processed. The Stent is placed on two moveable roles, due to the precise mechanics, the Stent can be placed precisely in position. The Stent Roller can be used with back and front light.



Technical Data

Maximal Length: 60 mm, 125 mm or 205 mm

Stent Ø: 0.5 mm, 1 mm or 6 mm

Standard Versions:

Articel Nr.	Туре	minimum Ø	maximum Length:	
27750	60-0,5	from 0,5 mm,	50 mm	
24130	125-1	from 1 mm,	125 mm	
27620	205-1	from 1 mm,	205 mm	
27590	125-6	from 6 mm,	125 mm	
27610	205-6	from 6 mm,	205 mm	

Miscellaneous:

Special sizes and customizations on request.



Made in Germany

All specifications are subject to change without notice.

Heinz Schade GmbH; Schießwieslenstraße 18; D-72766 Reutlingen